

ENGLISH

Indication for Use

Tendonitis, Carpal Tunnel Syndrome, arthritis, strains and sprains, post cast healing, and other selected injuries.

Caution

This product is designed to be worn directly on the skin. Inspect your hand, wrist, and forearm before and after each use. If any problems occur related to use of this product, please contact your healthcare provider.

Warning

Do not over-tighten the brace, risking loss of blood flow. Patients with already compromised blood circulation in their extremities, like those with diabetes or peripheral vascular disease need to be extra careful and are advised to consult with their doctor before using the brace. If you experience pain, swelling, sensation changes, or if your extremity shows signs of insufficient blood flow (turns blue, white or cold) while using this brace, discontinue the use of it immediately and consult with your medical professional. Never wear the brace directly over an open wound.